

Extraction of Dioxins and PCBs from Food and Feed

**Sean Panton, Frankie Smith, Mel Holland
Fera Science Ltd**

Principle

- The subsample is fortified with appropriate $^{13}\text{C}_{12}$ -labelled internal standards and subjected to exhaustive solvent extraction.
- The crude extract is mixed with solvent and sulphuric acid impregnated silica gel to remove co-extracted materials, including lipid and labile chemical contaminants.
- The slurry produced by this is quantitatively transferred to the large column upon which further extraction occurs.
- The eluent from the large column is passed through a column containing activated carbon. This retains the planar compounds of interest (PCDD/Fs, PBDD/Fs and non-ortho-PCBs/PBBs) whilst allowing the ortho-PCBs/PBBs and PBDEs to pass through.
- The planar compounds are eluted as a separate fraction



Method Schematic

Sample Preparation

Internal Standardisation of aliquot with ^{13}C -labelled dioxins and PCBs

50g / 5 g fat

Extraction and Fractionation

Hexane/Dichloromethane extraction/elution from Multilayered modified silica column

Fractionation on activated carbon

Other halogenated aromatics – PCBs PBDEs etc

Dioxins cPCBs

Purification

Dioxin fraction – purified on activated alumina

Concentration and Sensitivity Standard addition


Measurement

Analysis by HRGC-HRMS or GC-MS/MS

25 μl



Large Multilayered Column



Reagent Preparation

NaSO₄

- Washed with dichloromethane (2 x 3 L, then 4 L)
- Elute under N₂
- 120°C to 140°C for at least 12 hours
- Place in crystallising dishes in muffle oven for 5 to 9 days to remove water



Reagent Preparation

Silica Gel, activated

- Washed with dichloromethane (12.5 L then 10 L)
- Elute under N₂
- 130°C to 140°C for at least 12 hours
- Place in crystallising dishes in muffle oven for 5 to 9 days to remove water



Reagent Preparation

Silica Gel, Base modified

- KOH dissolved in methanol (280 g/L)
- Stir for 20 minutes
- Add silica gel
- Exothermic!!!
- Add methanol, mix slowly for 2 hours
- Elute methanol under N₂
- Elute with dichloromethane under N₂ (2 x 300 mL)



Reagent Preparation

Silica Gel, Acid modified

- Add 1 Kg H_2SO_4 to 1.5 kg activated silica
- Careful!!!
- Roller mix for at least 4 hours
- Free flowing powder



Reagent Preparation

Glass wool, Silanised

- Trimethylchlorosilane (5% v/v) in hexane
- Leave for 30 to 60 minutes
- Rinse with methanol then dichloromethane
- dry at 45°C to 50°C for at least 1 hour



Sample extraction

- Typically, 5g of fat or equivalent sample weight
- Reagent blank: 10 mL hexane
- Internal standard added
- Samples mixed with acid silica and hexane to create a slurry



Sample extraction

- Sample slurry added to top of multilayered column
- PCBs (and PBDEs) eluted with a 600 mL hexane: dichloromethane (60:40)
- Dioxins and non-ortho PCBs stick to carbon column



Sample Clean-Up

PCBs (and PBDEs)

- Only 1/5 of 600 mL extract used – concentrate to 1 mL
- Mini-column with NaSO₄ and acid silica
- Elute with 10 mL hexane
- Concentrate to 100 µL, add syringe standard reduce to 50 µL for analysis



Sample Clean-Up

Dioxins

- Reverse carbon column
- Elute with toluene (170 mL)
- Concentrate and exchange to hexane
- Further clean-up on activated alumina



Sample Clean-Up

Dioxins – activated alumina

- Elute with 40 mL hexane:dichloromethane (70:30)
- Concentrate to 100 μ L, add syringe standard
- reduce to 50 μ L for analysis



www.euchinasafe.eu



This project has received funding from the European Union's Horizon 2020 research and innovation programme under grant agreement No. 727864 and from the Chinese Ministry of Science and Technology (MOST).

Disclaimer: The content of this presentation does not reflect the official opinion of the European Commission and/or the Chinese government. Responsibility for the information and views expressed therein lies entirely with the author(s).

